

DenQ:

COMPANY
INTRODUCTION
REAFLET

Ver.2026



DenQ:
IMPLANT


**FDA 510k
Approved**



**MFDS
Approved**



Prosthetic Flow Diagram



ABOUT

DenQ is a company that develops a system of research and manufactures dental implants. Over 20 years of experience in manufacturing high-precision parts, and improving product quality. This system provides a wide range of products and solutions

COMPANY HISTORY

- 2019. 06. Establishment of DenQ
- 2019. 11. A Patent Application 'Dental Implant'
- 2020. 05. Registered 2 trademarks
- 2021. 07. A Patent Registration 'Dental Implant'
- 2021. 08. Obtained ISO 9001, ISO 13485 Certification
- 2022. 11. PCT International Patent Application 'Dental Implant Surface Treatment method'
- 2023. 10. A Patent Application 'Low-temperature Plasma surface treatment tool for CNC machine tools'
- 2023. 10. A Patent Application 'Cold Atmospheric Pressure Surface Modification Machine and Surface Modification Method using the same'
- 2024. 06. A Patent Registration 'Dental Implant Surface Treatment method'
- 2024. 09. Converted into a corporation
- 2024. 09. ISO 13485 Re-issue
- 2024. 09. FDA 510K Approved
- 2025. 08. MFDS Approved

Dental Implant System



DenQ Sub SLA Fixture

Economical Design
D3.7 through D5.1 All Hex 2.5 Regular
Connection

Typically emphasized features

Switching Platform
Double thread
SLA Surface Treatment

DenQ Taper KIT



CAD/CAM Components



ANALOG



Ti-BASE



SCAN BODY

Cement Retained Restoration



Cover Screw



Healing
Abutment



Temporary
Abutment



Cement
Abutment



Angled
Abutment



CCM Casting
Abutment



Fixture Level
Lab Analog



Pick-up
Impression Coping



Transfer
Impression Coping

Screw Retained Restoration



Healing Cap



Straight
Abutment



Angled
Abutment



Temporary
Cylinder



Ti Cylinder



CCM
Casting Cylinder



Lab Analog



Pick-up
Impression Coping




Transfer
Impression Coping

ENDO-REMOVER(G.P CUTTER) PATENT NO. 0317420

REMOVE EXCESS G.P IN THE BLINK OF AN EYE.



- EFFORTLESSLY REMOVE EXCESS G.P
- WITH THE ENDO-REMOVER



THIN AND LIGHTWEIGH.

The streamlined, lightweight design enhances comfort for both clinician and patient throughout the entire procedure.

SIMPLE AND SAFE.

Designed for user-friendly, a single-button interface allows the tip to reach the optimum temperature instantaneously. Maintenance is equally seamless, requiring only tip and battery replacement.

FOR LONG USE.

ENDO series are powered with rechargeable Ni-MH batteries, which are known to have 400~500 charge cycles on average.


The unit is compatible with standard replacement batteries available directly from us or via third-party providers. We offer a 1-year limited warranty from the date of purchase

CAUTION.

1. Before use, inspect the tip for any signs of damage.
2. Ensure correct battery by aligning the terminals with the (+) and (-) orientation markers engraved on the unit.
3. Ensure the batteries are charged before use.


HOW TO USE.

1. The tip heats up only when you press the button. Place the tip on the canal orifice.
2. Press the button for 2 seconds, and the tip will reach optimum temperature.
Remove excess G.P when the tip is heated.
3. Once you're finished, carefully extract the tip from the patient's mouth.




ENDO-CONDENSER UTILITY DESIGN APPLICATION

↙ HIGH-DENSITY G.P CHARGE, PERFECT OBTURATION.



- EASILY CHARGE GUTTA-PERCHA
- OBTURATE ROOT CANAL WITH ENDO-CONDENSER



THIN AND LIGHTWEIGH.

The ENDO-CONDENSER is only 1cm thick and weighs just 58 grams, batteries included. It's designed to be lightweight, making the entire treatment process more comfortable.

PRECISE CHARGING & CONDENSATION, SUCCESSFUL TREATMENT.

Press the button for 2 seconds to achieve the perfect temperature for softening Gutta-Percha. It makes filling and obturating the root canal simple, fast, and accurate. Get a denser, more even G.P. compaction every time. The ENDO-CONDENSER is a true multi-tool: use it as a spreader, plugger, or heat carrier for all your warm vertical condensation needs.

FOR LONG USE.

Maintenance is easy; you may replace the batteries and tips as required. Moreover, this unit comes with a 1-year warranty and free repairs from the date of purchase.

CAUTION.

1. Please inspect the tip for any damage before use to ensure the best results
2. Make sure to insert the batteries correctly, following the + and - signs printed on the unit.
3. Ensure the batteries are charged before each use.

HOW TO USE.

1. The tip heats up only when you press the button. Before treatment, you may warm up the device for about 5 seconds. Once ready, insert the tip and press the button to fill the G.P. material into your desired position.
2. Make sure to hold the button for about 4–5 seconds before you remove the tip from the canal.
3. Once the treatment is finished, let the tip cool down completely before sanitizing it.



ORDER INFORMATION



ENDO-REMOVER

REMOVER 1EA
CHARGER 1EA
INSTRUCTIONS FOR
USE AND WARRANTY

BATTERY CELL 2EA
CUTTING TIPS 4EA



ENDO-CONDENSER

CONDENSER 1EA
CHARGER 1EA
INSTRUCTIONS FOR
USE AND WARRANTY


BATTERY CELL 2EA
TIP 2EA : FINE,
FINE MEDIUM

BITE(EXCEL) TRAY

RELIABLE IMPRESSION TAKING.



- Rigid Frame
- Ultra-Thin Contact
- Thin and Durable Non-woven Fabric
- Easy, Precise Impression



ORDER INFORMATION

- A. POSTERIOR : 50 Trays + 50 Stickers PER RACK
- B. SIDELESS POSTERIOR : 50 Trays + 50 Stickers PER PACK
- C. FULL ARCH : 25 Trays + 25 Stickers PER PACK
- D. ANTERIOR : 30 Trays + 30 Stickers PER PACK
- E. QUADRANT : 30 Trays + 30 Stickers PER PACK

<5 SIZES> STICKERS (DATE, NAME, AGE, SEX) INCLUDED.
ABUTMENT, ANTAGONIST TEETH, DENTAL PROSTHESIS
IMPRESSION. AT THE SAME TIME.

MIXING TIP

EASIER MIXING AND INJECTION.



HEAVY



LIGHT

MIXTIP MADE OUT
OF RANDOM PP

ORAL TIP MADE OUT
OF RANDOM PP

NOW EASILY MIX & INJECT IMPRESSION MATERIAL WITH MIXING TIPS.
MEET THE HIGH-QUALITY DENTAL SEARCH MIXING TIPS AT A REASONABLE PRICE.


*The color of the product on each photo depends on the lighting, camera settings, etc.
The actual color is most similar to the photo below

ORDER INFORMATION

- HEAVY : 48 Tips PER RACK
- LIGHT : 48 Tips & 20 Oral Tips PER PACK
- Oral Tips Only : 100 PER PACK



EXCEL PLASTIC TRAY FOR BETTER IMPRESSION TECHNIQUE



RIGID FRAME STRUCTURE.

High-strength PS Frame.

OPTIMAL WALL HEIGHT.

Reduced Pressure on the Gums.

Maximizing Operator and Patient Comfort.

10 SIZE OPTIONS

or various clinical applications.

MINIMUM DISCOMFORT

in the oral cavity.

ORDER INFO

POSTERIOR



PB(#1-1)
48EA / PACK

PA(#1-2)
48EA / PACK

ANTERIOR



AU(#2-1)
48EA / PACK

AL(#2-2)
48EA / PACK

FULL ARCH



UL(#3-1)
40EA / PACK

LL(#3-2)
40EA / PACK

UM(#4-1)
48EA / PACK

LM(#4-2)
48EA / PACK

US(#5-1)
48EA / PACK

LS(#5-2)
48EA / PACK



DenQ Implant
WhatsApp Biz

INFORMATION

+82 10 8210 9792

biz@denq.kr www.denq.kr

302 Hyowon, Academic Industry 2, Busandaehak-ro 63beon-gil, Geumjeong-gu, Busan, Korea